



SUITABILITY TESTING BETWEEN COMMERCIAL SPE SORBENTS FOR THE SAMPLE CLEAN-UP IN STA USING LC-MS(/MS)

Tineke Decaestecker

E. Coopman, C. Van Peteghem, and J. Van Bocxlaer

Laboratory of Medical Biochemistry and Clinical Analysis,
Ghent University, Belgium



Suitability testing between commercial SPE sorbents for the sample clean-up in STA using LC-MS(/MS)



Outline

- Introduction and goal
 - ✓ DDA - fundamentals
- Study set-up
 - ✓ Flowchart
 - ✓ Extraction yield evaluation
 - ✓ Clean-up potential evaluation
- Results
 - ✓ Evaluated SPE-sorbents
 - ✓ Extraction yield evaluation
 - ✓ Clean-up potential evaluation
- Conclusion



Introduction & Goal

- **General Unknown Screening procedures**

 - ☞ no foreknowledge

 - ⇒ **Data-Dependent Acquisition (DDA) ~ Q-TOF***

 - * RCM 14, p. 1787-1792 (2000)

- **Suitability test of SPE sorbents:**

 - ✓ Apolar

 - ✓ Mixed-mode

 - ✓ Polymeric

extraction yield ↔ cleanliness




Selection of best sorbent


DDA – fundamentals (1)

- Initially: QUAD = wide band-pass filter
 - ◆ precursor ion(s) > MS threshold
 - ⇒ switch to MS/MS
 - fragment ions ⇒ TOF
 - ◆ fragment ion(s) < MS/MS threshold
or total MSMS acquisition time
 - ⇒ switch back to MS

DDA – fundamentals (2)




DDA – fundamentals (3)



Advantages DDA

- No foreknowledge
 - Quantitative indication
 - Clean MS/MS spectra
- } 1 run

Study set-up: flowchart



Extraction yield evaluation

For sorbent_x

	EY	
Compound 1	A%	Value 1 = F(A%)
Compound 2	B%	Value 2 = F(B%)
.	.	.
.	..	.
Compound 18	R%	Value 18 = F(R%)
		<hr/>
		V = \sum Value 1..18


F() is a particular transformation function (see next slide).

Comparing the extraction yield of the sorbents

=

comparing the value V of the sorbents

Transformation function for extraction yield evaluation



Characteristics:

if **Yield < 50%**: more than linear **penalty**


if **Yield > 50%**: more than linear **reward**

Clean-up potential evaluation

- 2 criteria:
 - ✓ DDA: number of spiked compounds “recovered”
 - ✓ CV% of EY

Results: 12 sorbents (3 categories)

APOLAR



MIXED-MODE

HCX: C8/SO₃⁻

HCX3: C18/ SO₃⁻

HCX5: C4/ SO₃⁻

POLYMERIC

HLB: DVB-NP


MCX: DVB-NP/ SO₃⁻

R=(CH₂)_nCH₃ C2, C4, C8, C18, C18MF

R=C₆H₅ PH


R=(CH₂)₃CN CN

Results: extraction yield evaluation




- Evaluation of 12 sorbents based on several extractions in triplicate

Results: extraction yield evaluation



- Evaluation of 12 sorbents based on several extractions in triplicate
- Choice of 2 sorbents for clean-up potential evaluation: C8 and MCX

Results: clean-up potential evaluation



- Fragmentation of 15 compounds
- Choice of sorbent: C8

Conclusions

- Evaluation of a range of SPE-sorbents on extraction yield and cleanliness criteria
- Selection of C8 for future use with DDA
- DDA: useful for STA